# **Dynamixel Wizard Configuration**

- 1) Download the package file.
  - a. Linux Software Download
- 2) Enter the following command to change the permission.

### \$ sudo chmod 775 DynamixelWizard2Setup\_x64

3) Run the install program.

# \$ ./DynamixelWizard2Setup\_x64

- 4) Click on Next button to proceed installation.
- 5) After completing the installation, please add account id to dialout group in order to access the USB port. Replace the in the command below to your actual user id.

## \$ sudo usermod -aG dialout <your\_account\_id>

6) Reboot for the changes to be effective.

### \$ reboot

#### Once it is installed

- 1. Scan the ports for the dynamixels
- 2. Once connected, Turn on the **Torque** for each motor and play around with the angle knob to see the robot reaches different positions
- 3. Remember to close the Wizard once you start working with ROS2 Control

#### Tips:

- The USB root permissions are important as the motors will not be recognized during scanning otherwise
- 2. Most times, a power reset/ unplugging can solve the issue

### **Open Manipulator X with ROS2**

This part with ROS2 has a prerequisite of creating and building your workspace with the following packages provided as part of the assignment.

- 1. open\_manipulator\_x
  - a. open\_manipulator
  - b. open\_manipulator\_x\_controller
  - c. open\_manipulator\_x\_description
  - d. open\_manipulator\_x\_libs
  - e. open\_manipulator\_x\_dependencies
  - f. openmanip\_control
  - g. open\_manipulator\_msgs
  - h. robotis\_manipulator
- 2. rbe500-example, dynamixel\_sdk\_examples *Optional but recommended for the assignment*

Once the workspace is built, source your workspace and run the following command to check if the controller runs and connects to the robot

\$ ros2 launch open\_manipulator\_x\_controller open\_manipulator\_x\_controller.launch.py

Now, check the ros2 topic list and you should be able to see /joint\_states

Echo the /joint\_states topic and you will find the joint angles of each of the 4 joints and the gripper joint in realtime.

You can also visualize the robot in rviz2 with

\$ ros2 launch open\_manipulator\_x\_description open\_manipulator\_x\_rviz.launch.py

The services used to control the robot can be found <u>here</u>. These are the services you will call to make your joints move to desired positions.

You can run this command to call a service to set the joint angles:

\$ ros2 service call /goal\_joint\_space\_path open\_manipulator\_msgs/srv/SetJointPosition "{joint\_position: {joint\_name: ['joint1', 'joint2', 'joint3', 'joint4'], position: [1.0, -1.0, 0.5, 0.0], max\_accelerations\_scaling\_factor: 0.5, max\_velocity\_scaling\_factor: 0.5}}"

For more information and some troubleshooting help, refer to the Official Robotis page